



## APPLICATIONS MANUAL

### CONTENTS

Ifosfamide in blood serum .....	2
Ingredients in blood serum .....	3
Organophosphorus pesticides in tea leaf .....	5
Sudan i ii iii iv in chilli sauce .....	6
Pah in water .....	7
Phenols in water .....	8
Tricyclic antidepressants in serum and plasma .....	9
Clonazepam & 7-aminoclonazepam in urine .....	10
Gabapentin in serum, plasma, or whole blood .....	11
Amphetamines in urine .....	12
Cocaine/benzoylecgonine in urine .....	13
Opiates in urine .....	14
Carboxy-thc in urine .....	15



## IFOSFAMIDE IN BLOOD SERUM

### CHROMPURE™ C18 SPE COLUMN (Cat# LBSC181001)

#### PREPARE SAMPLE

Solution 0.2g ifosfamide samples in 100ml volumetric flask with mobile phases to be 2mg/ml internal standard, add 12.5ul internal standard in 0.5ml serum samples.

#### CONDITION COLUMN

condition SPE column with 2ml Acetonitrile and 2mL Physiological saline.

#### PURIFICATION


load the serum samples into SPE columns (chrompure C18)

Wash the SPE column with 1ml Physiological saline and followed by washing with 5% Acetonitrile, concentrate to dryness with N<sub>2</sub>

Elute with 0.5ml Acetonitrile

#### CHROMATOGRAPHY

Inject 20 µl onto UV-HPLC, flow rate 1.0ml/min, room temp , det. UV-Vis 200nm



Concentration(ug/ml)	Recoveries(%)
100	92.0
50	93.5
5	89.3



## INGREDIENTS IN BLOOD SERUM

### CHROMPURE™ PLS SPE COLUMN (Cat# LBSPLS5003)

#### PREPARE SAMPLE

To 1 mL of serum or plasma add internal standard\* and 2 ml of 100 mM phosphate buffer (pH 6.0). Mix/vortex.

Centrifuge for 10 minutes at 2000 rpm and discard pellet. Sample pH should be  $6.0 \pm 0.5$ . Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.

#### CONDITION COLUMN

condition SPE column with 3ml Methanol and 3ml H<sub>2</sub>O

#### PURIFICATION



load the serum samples into SPE columns (chrompure PLS)

Wash the SPE column with 3ml H<sub>2</sub>O or 5% Methanol

Elute with 3ml Methanol

#### CHROMATOGRAPHY

Inject 20 µl onto UV-HPLC, flow rate 1.0ml/min, room temp .



Matter	Recoveries(%)
Dexamethasone	97.9
Ethinylestradiol	96.3
Hydrocortison $\mu\text{m}$	74.0
Triamcinolone	71.9
Levonorgestrel	93.9
Ganciclovir	54.1
Prednisone Acetate	98.6
Cefalexin	58.6
Cefradine	45.6



## **ORGANOPHOSPHORUS PESTICIDES IN TEA LEAF CHROMPURE™ GCB SPE COLUMN (CAT#LBSGCB5006)**

### **PREPARE SAMPLE**

2g tea leaf add into water over night. Homogenize the tea leaf in 10ml(n-Hexan:acetone=1:1), repeat homogenize for 2-3 times. Unite all the liquid, homogenize, centrifuge and get the upper layer, Add proper anhydrous sodium sulfate and centrifuge , get the upper layer, evaporate supernatant to 1 mL under N<sub>2</sub> at 35 °C

### **CONDITION COLUMN**

5 mL n-Hexan:acetone=1:1, GCB 500mg/6ml (Cat#LBSGCB5006)

### **PURIFICATION**

Load 1ml sample onto GCB SPE column

elute with 20ml n-Hexan:acetone=1:1

Gather all samples and concentrate to dryness with N<sub>2</sub> at 40 °C,

Reconstitute in 1 mL acetone:n-hexane (1:1)

### **CHROMATOGRAPHY**

Inject 1.0 µl onto GC-FPD 50 °C (1min), 15 °C /min to 200 °C (2min), 5 °C /min to 260 °C (8min), FPD,250 °C



## SUDAN I II III IV IN CHILLI SAUCE

**CHROMPURE™ Alumina-N SPE COLUMN (Cat# LBSALNB5003)**

### PREPARE SAMPLE

Homogenize 5g chilli sauce in 10ml(n-Hexan:acetone=3:1), ultrasonic 15min, centrifuge , get the acetone layer, ultrasonic the residual with 5ml n-Hexan for 2 time, gather the acetone layer, Add proper anhydrous sodium sulfate, filter by 0.45um membrane,concentrate to dryness with N<sub>2</sub> to 5ml

### CONDITION COLUMN

5 mL n-Hexan, Al-N 500mg/3ml(Cat# LBSALNB5003)

### PURIFICATION

Load the sample onto Al-N SPE column at speed 1 second/drop

Wash with 3\*5ml n-Hexan to dry

Elute with 5ml n-Hexan(5% acetone)

Concentrate to dryness with N<sub>2</sub> at 40 °C, fix to 5ml by methanol

### CHROMATOGRAPHY

Inject 20 µl onto UV-HPLC, flow rate1.0ml/min, temp 35 °C, det. UV-Vis  
500nm



## PAH IN WATER

### CHROMPURE™ C18 SPE COLUMN (Cat# LBSC1810006)

#### PREPARE SAMPLE

6N HCL add to 250ml – 1L water adjust to PH<2

#### CONDITION COLUMN

6 mL dichloromethane, 6ml methanol, 6ml deionized water, C18 1000mg/6ml (Cat# LBSC1810006)

#### PURIFICATION

Load the water on into SPE tube 5ml/minute

Wash C18 column with 6ml deionized water, dry under vacuum pressure

Elute by 3\*1ml dichloromethane

Gather all sample, concentrate to dryness with N<sub>2</sub> at room temp to 0.5ml, add perylene-d12 to be internal standard

#### CHROMATOGRAPHY

Inject 1.0 µl onto GC-MSD, at 100 °C (1min), 6 °C /min to 300 °C (30min), MSD280 °C, carrier gas: He 1.0ml/min, m/z 50-450



## **PHENOLS IN WATER**

### **CHROMPURE™ PLS SPE COLUMN (Cat# LBSPLS1003)**

#### **PREPARE SAMPLE**

100mL water spiked with phenols, if sample containing in sediment may require prefiltration.

#### **CONDITION COLUMN**

6 mL methyl t-butyl ether, 6 mL methanol, 6 mL deionized water, PLS 100mg/3ml(Cat# LBSPLS1003)

#### **PURIFICATION**

Load 100 mL aqueous sample onto SPE tube in less than 5ml/min

Elute with 2\*2.5mL methyl t-butyl ether and allow elution to proceed at a dropwise rate, add methyl t-butyl ether to 5ml

#### **CHROMATOGRAPHY**

Inject 1.0 µl onto GC-MSD, 200 °C, splitless (45 sec delay), temp. at 65 °C to 185 °C at 10 °C/min, hold 1 min, then to 275 °C, at 20 °C/min, hold 5 min. FID, 330 °C, Carrier gas: Nitrogen, 1.0 mL/min





## **TRICYCLIC ANTIDEPRESSANTS IN SERUM AND PLASMA CHROMPURE™ PXA SPE COLUMN (Cat# LBSPAX1003)**

### **PREPARE SAMPLE**

To 1 mL of serum or plasma add internal standard\* and 2 ml of 100 mM phosphate buffer (pH 6.0). Mix/vortex. Centrifuge for 10 minutes at 2000 rpm and discard pellet. Sample pH should be  $6.0 \pm 0.5$ . Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.

### **CONDITION COLUMN**

3 mL CH<sub>3</sub>OH. 3 mL D.I. H<sub>2</sub>O. 1 mL 100 mM phosphate buffer (pH 6.0). Aspirate at < 3 inches Hg to prevent sorbent drying. PAX 200mg/3ml(Cat# LBSPAX2003)

### **PURIFICATION**

Load at 1 mL/minute. Wash with 3 mL D.I. H<sub>2</sub>O. 1 mL 100 mM acetic acid. 3 mL CH<sub>3</sub>OH. Dry column (5 minutes at > 10 inches Hg)

Elute with 3 mL CH<sub>2</sub>Cl<sub>2</sub>/IPA/NH<sub>4</sub>OH (78:20:2). Collect eluate at 1 mL/minute or use gravity flow. Evaporate to dryness at < 40 °C.

### **CHROMATOGRAPHY**

Inject 100 µL onto HPLC. Column temp = 30 °C Mobile phase- Acetonitrile/ Buffer/ Methanol (60:25:15), Buffer = 0.01 M K<sub>2</sub>HPO<sub>4</sub> adjusted to pH 7.0 with H<sub>3</sub>PO<sub>4</sub> Flow rate = 1.75 mL/min.



## **CLONAZEPAM & 7-AMINOCLONAZEPAM IN URINE CHROMPURE™ PLS SPE COLUMN (Cat# LBSPLS0301)**

### **PREPARE SAMPLE**

To 2 mL of urine, add internal standard(s)\* and 1 mL of  $\beta$ -Glucuronidase solution.  $\beta$ -Glucuronidase solution contains 5,000 F units/mL *Patella vulgata* in 100 mM acetate buffer (pH 5.0). Mix/vortex. Hydrolyze for 3 hours at 65 °C. Cool before proceeding.

### **CONDITION COLUMN**

3 mL CH<sub>3</sub>OH. 3 mL deionized water. 1 mL 100 mM phosphate buffer (pH 6.0).

### **PURIFICATION**

Load samples at 1 to 2 mL/ minute.

Wash with 2 mL deionized water. 2 mL 20% acetonitrile in 100 mM phosphate buffer (pH 6.0). Dry column (5 minutes at > 10 inches Hg). 2 mL hexane.

Elute with 3 mL ethyl acetate with 2% NH<sub>4</sub>OH: Collect eluate at 1 to 2 mL/minute.

Evaporate to dryness at < 40 °C.

Add 50  $\mu$ L ethyl acetate and 50  $\mu$ L MTBSTFA(with 1% TBDMCS), Mix/vortex.

React 20 minutes at 90 °C. Remove from heat source to cool.

### **CHROMATOGRAPHY**

Inject 1 to 2  $\mu$ L sample



## GABAPENTIN IN SERUM, PLASMA, OR WHOLE BLOOD CHROMPURE™ C18 SPE COLUMN (Cat# LBSC181001)

### PREPARE SAMPLE

500 µL of sample, calibrator, or control was placed into a disposable glass test tube and 25 µL of internal standard\* (5.0 mg/L) was added. Vortex tube. Add 500 µL of 20% acetic acid and vortex tube again.

### CONDITION SPE COLUMN

1 x 3 mL CH<sub>3</sub>OH.  
1 x 3 mL D.I. H<sub>2</sub>O.  
1 x 1 mL 100 mM HCL.

### APPLY SAMPLE

Load at 1 to mL/minute.

### WASH COLUMN

1 x 3 mL D.I. H<sub>2</sub>O.  
1 x 3 mL ethyl acetate.  
1 x 3 mL hexane.  
Dry column.  
(5 minutes at > 10 inches Hg) or until column is dry.

### ELUTION

1 x 1 mL 2% NH<sub>4</sub>OH in CH<sub>3</sub>OH.

### DRY ELUATE

Evaporate to dryness at < 40 °C

### DERIVATIZATION

Add 50 µL of MTBSTFA + 1 % t-BDMCS\*\* and 50 µL ethyl acetate.  
Cap and heat at 70 °C for 30 minutes.  
Remove and allow to cool.

### QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.



## AMPHETAMINES IN URINE CHROMPURE™ PCX SPE COLUMN (Cat# LBSPCX301)

### SAMPLE PREPERATION

To 1 mL of urine add internal standard(s) and 1 mL 100mM phosphate buffer (pH 6.0). Mix/Vortex.

### APPLY SAMPLE TO COLUMN

Load at a rate of 1 to 2 mL/min.

### WASH COLUMN

1 x 1 mL DI H<sub>2</sub>O.

1 x 1 mL 100mM acetic acid.

1 x 1 mL MeOH.

Dry column (3 mins at > 10 inches Hg).

### ELUTE AMPHETAMINES

2 x 0.5 mL CH<sub>2</sub>Cl<sub>2</sub>/IPA/NH<sub>4</sub>OH (78/20/2), collect eluate at 1 to 2 mL/min.

### CONCENTRATE ELUATE

Add 1 drop 1% HCl in MeOH to eluate before evaporating.

Evaporate to dryness at < 40 °C.

### DERIVATIZATION

Add 50 ul ethyl acetate and 50 ul TFA (Trifluoroacetic acid anhydride) then cap, mix/vortex.

Heat for 15 mins at 70 °C, allow to cool, then evaporate to dryness at < 40 °C.

Reconstitute with 100 µL ethyl acetate.

### ANALYZE

Inject 1 to 2 µL onto gas chromatograph. For MSD monitor the following ions:

Analyte (TFA)	Target (Quantitation) Ion	Qualifier Ions
Amphetamine	140	91, 118
Amphetamine-d11*	144	98, 128
Methamphetamine	154	110, 118
Mehtamphetamine-d11*	160	113, 126

\*Suggested internal standards



## COCAINE/BENZOYLECGONINE IN URINE CHROMPURE™ PCX SPE COLUMN (Cat# LBSPCX301)

### SAMPLE PREPERATION

To 1 mL of urine add internal standard(s) and 300 µl 100mM HCl.  
Mix/Vortex.

### APPLY SAMPLE TO COLUMN

Load at a rate of 1 to 2 mL/min.

### WASH COLUMN

1 x 1 mL DI H<sub>2</sub>O.

1 x 1 mL 100mM HCl.

1 x 1 mL MeOH.

Dry column (3 mins at > 10 inches Hg).

### ELUTE COCAINE/BENZOYLECGONINE

2 x 0.5 mL CH<sub>2</sub>Cl<sub>2</sub>/IPA/NH<sub>4</sub>OH (78/20/2), Collect eluate at 1 to 2 mL/min.

### CONCENTRATE ELUATE

Evaporate to dryness at < 40 °C.

### DERIVATIZATION

Add 50 µL ethyl acetate and 50 µL BSTFA w/ 1% TMCS, then cap, mix/vortex. Heat for 20 mins at 70 °C, allow to cool.

### ANALYZE

Inject 1 to 2 µL onto gas chromatograph. For MSD monitor the following ions:

Analyte	Target (Quantitation) Ion	Qualifier Ions
Cocaine	182	198, 303
Cocaine-d3*	185	201, 306
Benzoylecgonine (TMS)	240	256, 361
Benzoylecgonine-d8 (TMS)*	243	259, 369

\*Suggested internal standards



## OPIATES IN URINE CHROMPURE™ PCX SPE COLUMN (Cat# LBSPCX301)

### SAMPLE PREPERATION

To 1 ml of urine add internal standard(s) and 1.0 ml  $\beta$ -Glucuronidase solution. ( $\beta$ -Glucuronidase solution contains 5000 Funits/mL *Patella Vulgata* in 100mM acetate buffer, pH 5.0). Hydrolyze for 3 hours at 60 °C. Cool, then centrifuge for 10 minutes at high speed and discard pellet. Adjust pH to 6.0  $\pm$ 0.5 with 1.0N NaOH. NOTE: For unconjugated (free) opiates; to 1 mL urine, add internal standard(s) and 1 mL 100mM phosphate buffer (pH 6.0). Proceed to Step #2.

### APPLY SAMPLE TO COLUMN

Load at a rate of 1 to 2 mL/min.

### WASH COLUMN

- 1 x 1 mL DI H<sub>2</sub>O.
- 1 x 1 mL 100mM acetate buffer (pH 4.5).
- 1 x 1 mL MeOH.
- Dry column (3 mins at > 10 inches Hg).

### ELUTE OPIATES

2 x 0.5 mL CH<sub>2</sub>Cl<sub>2</sub>/IPA/NH<sub>4</sub>OH (78/20/2), collect eluate at 1 to 2 mL/min. Evaporate eluate to dryness at < 40 °C.

### DERIVATIZATION

Add 50  $\mu$ L ethyl acetate and 50  $\mu$ L BSTFA w/ 1% TMCS, then cap, mix/vortex. Heat for 20 mins at 70 °C, allow to cool.

### ANALYZE

Inject 1 to 2  $\mu$ L onto gas chromatograph: For MSD monitor the following ions:

Analyte (TMS)	Target (Quantitation) Ion	Qualifier Ions
Codeine	371	234, 343
Codeine-d6*	377	237, 349
Morphine	429	401,414
Morphine-d6*	435	404, 420
6-Acetylmorphine	399	287, 340

\*Suggested internal standards



## CARBOXY-THC IN URINE CHROMPURE™ PCX SPE COLUMN (Cat# LBSPCX301)

### SAMPLE PREPERATION

To 2 mL of urine add internal standard and 100 µL 10N NaOH. Mix/vortex. Hydrolyze for 20 mins at 60 °C. Cool before proceeding. Adjust sample pH to 3.5 ±0.5 with 1.0 mL glacial acetic acid.

### APPLY SAMPLE TO COLUMN

Load at a rate of 1 to 2 mL/min.

### WASH COLUMN

1 x 1 mL DI H<sub>2</sub>O.  
1 x 1 mL 0.1M HCl/acetonitrile (70/30).  
Dry column (3 mins at > 10 inches Hg).  
1 x 200 µL hexane.

### ELUTE CARBOXY-THC

2 x 0.5 mL hexane/ethyl acetate (50:50); Collect eluate at 1 to 2 mL/min.  
Evaporate eluate to dryness at < 40 °C.

### DERIVATIZATION

Add 50 µL ethyl acetate and 50 µL BSTFA w/ 1% TMCS, then cap, mix/vortex.  
Heat for 20 mins at 70 °C, allow to cool.

### ANALYZE

Inject 1 to 2 µL onto gas chromatograph. For MSD monitor the following ions:

Analyte (TMS)	Target (Quantitation) Ion	Qualifier Ions
Carboxy-THC	371	473, 488
Carboxy-THC-d3*	374	476, 491

\*Suggested internal standards