




## MS® Serological Pipette-GRADUATED

*ISO9001: 2000; FDA registered; CE certification*



MS® Serological Pipette fulfil all requirements set for disposable serological pipettes: **crystal clear** transparency, bidirectional graduations that provide a choice of ascending and descending volume reading, and accurate anti-drip delivery.

The cotton wadding sits firmly in the tube hinders an over-filling of the pipette sucking apparatus. The **colour**-coding on the pipettes assists in avoiding confusion during routine examinations and make it simple to determine the size/volume of each pipette.

Pipettes are pyrogen free, non cytotoxic, non haemolytic and come single wrapped in a convenient blister package. All the features guarantee them can be safely used in cell and tissue culture, bacteriological, and clinical research.



Sterile Serological Pipette



The same from their back

Items	Capacity/ Subdivisions (ml)	Color Code	Sterile
LBGSP000001	1 / 0.01	YELLOW	Yes
LBGSP000002	2 / 0.01	GREEN	Yes
LBGSP000005	5 / 0.1	BLUE	Yes
LBGSP000010	10 / 0.1	RED	Yes
LBGSP000025	25 / 0.2	PURPLE	Yes
LBGSP000050	50 / 1.0	BLACK	Yes

Certified as having been sterilized by gamma radiation

Items	Capacity/ Subdivisions (ml)	Color Code	Sterile
LBGSP001001	1 / 0.01	YELLOW	NO
LBGSP001002	2 / 0.01	GREEN	NO
LBGSP001005	5 / 0.1	BLUE	NO
LBGSP001010	10 / 0.1	RED	NO
LBGSP001025	25 / 0.2	PURPLE	NO
LBGSP001050	50 / 1.0	BLACK	NO

**Membrane Solutions, LLC**