


MS[®] Serological Pipettes

Stay Pure



MS® Serological pipettes are the most favorable product of liquid handling for life science products profile. All the serological pipettes have been manufactured under 100K grade cleanroom cleanness with ISO 9001:2000 and CE certificated.

High Precision


- ◇ Exact graduation (<2%)
- ◇ Inside slick surface
- ◇ Drop-free tip design
- ◇ Strict leakage tested

Easy to use

- ◇ 100% virgin polystyrene, highly visible
- ◇ Customary standard length
- ◇ Cotton filter plug for saving Pipet-Aids
- ◇ Easy-to-rip paper-plastic pack

Sterility

- ◇ Non-pyrogenic
- ◇ Sterilized by gamma-radiation with independent packing



Order Information

| Item | Total volume (mL) | Graduation (mL) | Negative Grads (mL) | Color code | Qty/Pk | Qty/Cs |
|---------|-------------------|-----------------|---------------------|------------|--------|--------|
| LBSP01S | 1 | 1/100 | 0.4 | Yellow | 500 | 3000 |
| LBSP02S | 2 | 1/100 | 0.4 | Green | 400 | 2400 |
| LBSP05S | 5 | 1/10 | 2.5 | Blue | 300 | 1800 |
| LBSP10S | 10 | 1/10 | 3 | Orange | 200 | 1200 |
| LBSP25S | 25 | 2/10 | 8 | Red | 200 | 800 |
| LBSP50S | 50 | 1/2 | 10 | Black | 100 | 600 |

For more information, please contact:

USA

Fax: 732-412-4040

Address: 1108 Spring view Lane, Plano, TX 75075

China

Toll free: 400-820-0293

Fax: 0086-21-51687551

Address: 2202, No.1759 North Zhongshan Road, Shanghai

Japan

Toll free: 0066-33-800658 Toll free: 1-866-528-4572

Address: Suite #504, 1-19-2 Nihonbashi, kakigaracho, Chou-ku, Tokyo 103-0014

Inquiry

Please write to info@membrane-solutions.com

<http://www.membrane-solutions.com>



Membrane Solutions