


MS- Vials



2mL Screw Top Vials



LBSV082C	2mL clear vial, 10-425 screw top	100/Pk
LBSV002C	2mL clear vial, 9-425 screw top	100/Pk
LBSV012C	2mL clear vial, 8-425 screw top	100/Pk
LBSV092C	2mL clear vial, 10-425 screw top, graduated with writing area	100/Pk
LBSV022C	2mL clear vial, 9-425 screw top, graduated with writing area	100/Pk
LBSV032C	2mL clear vial, 8-425 screw top, graduated with writing area	100/Pk
LBSV082A	2mL amber vial, 10-425 screw top	100/Pk
LBSV002A	2mL amber vial, 9-425 screw top	100/Pk
LBSV012A	2mL amber vial, 8-425 screw top	100/Pk
LBSV092A	2mL amber vial, 10-425 screw top, graduated with writing area	100/Pk
LBSV022A	2mL amber vial, 9-425 screw top, graduated with writing area	100/Pk
LBSV032A	2mL amber vial, 8-425 screw top, graduated with writing area	100/Pk



2mL Crimp Top Vials



LBSV042C	2mL clear vial, crimp top	100/Pk
LBSV062C	2mL clear vial, crimp top, graduated with writing area	100/Pk
LBSV042C	2mL amber vial, crimp top	100/Pk
LBSV062C	2mL amber vial, crimp top, graduated with writing area	100/Pk

Septa for 2mL Vials




LBSV02RS	White PTFE/red Rubber septa, for 2mL 8-425 screw top vial	10000/Pk
LBSV02SS	Red PTFE/white silicone septa, for 2mL 8-425 screw top vial	10000/Pk
LBSV03RS	White PTFE/red Rubber septa, for 2mL 9-425 screw top vial	10000/Pk
LBSV03SS	Red PTFE/white silicone septa, for 2mL 9-425 screw top vial	10000/Pk
LBSV23SS	Blue PTFE/white silicone septa, Pre-slit, for 2mL 9-425 screw top vial	10000/Pk
LBSV33SS	Red PTFE/white silicone/red PTFE septa, for 2mL 9-425 screw top vial	10000/Pk
LBSV04RS	White PTFE/red Rubber septa, for 2mL crimp top vial	10000/Pk
LBSV04SS	Red PTFE/white silicone septa, for 2mL crimp top vial	10000/Pk



Cap for 2mL Vials


LBSV002	Blue screw cap with hole, for 2mL 9-425 screw top vial	100/Pk
LBSV012	Black screw cap with hole, for 2mL 8-425 screw top vial	100/Pk
LBSV032	Blue screw cap with hole, for 2mL 8-425 screw top vial	100/Pk
LBSV022	Aluminium cap with hole, for 2mL crimp top vial	100/Pk



Septa + Cap for 2mL Vials


LBSV002CRS	White PTFE/red Rubber septa + Blue screw cap with hole, for 2mL 9-425 screw top vial	100/Pk
LBSV012CRS	White PTFE/red Rubber septa + Black screw cap with hole, for 2mL 8-425 screw top vial	100/Pk
LBSV102CSS	Red PTFE/white silicone septa + Blue screw cap with hole, for 2mL 9-425 screw top vial	100/Pk
LBSV012CSS	Red PTFE/white silicone septa + Black screw cap with hole, for 2mL 8-425 screw top vial	100/Pk
LBSV112CRS	White PTFE/red rubber septa + Aluminium cap with hole, for 2mL crimp top vial	100/Pk
LBSV122CRS	Red PTFE/ white silicone septa + Aluminium cap with hole, for 2mL crimp top vial	100/Pk
LBSV302CRS	Blue PTFE/white Rubber septa, pre-slit + Blue screw cap with hole, for 2mL 9-425 screw top vial	100/Pk
LBSV402CSS	PTFE/ silicone/PTFE septa+ Blue screw cap with hole, for 2mL 9-425 screw top vial	100/Pk

2mL crimp top vials




d1	11.6±0.3
h1	32.0±0.5
d2	10.7±0.2
d4	6.0±0.2
h4	3.6±0.2
d3	8.3±0.3
h3	≈7.5
S2	≥0.4
S1	0.85±0.05

2mL screw top vials (8-425)



d1	11.6±0.25
d2	≤8.0
d4	5.0±0.30
h1	32.0±0.50
d3	8.8±0.30
h3	≈7.8
h5	1.6
h6	1.0
S2	≥0.50
S1	0.85±0.04

2mL screw top vials (9-425)



d1	11.6±0.2
d2	8.1-8.5
d4	6.0±0.20
h1	32.0±0.50
d3	9.4±0.25
h3	7.3-7.9
h5	≈1.2
h6	≈0.9
S2	≥0.50
S1	0.85±0.04

4mL, 10mL, 20mL, 30mL, 40mL sample vials and More Product

For More Information, Please email to: info@membrane-solutions.com